

InjectableRF® Recommended Supplies

Equipment and Disposable Items:

- Thermi® Radiofrequency Generator and cart
- FLIR® Camera, Monitor, Boom
- Nerve mapping device such as Xavant Stimpod
- Thermi InjectableRF® Electrode Disposable (*Size selected by physician for treatment area*)
 - AC-SN-DUS (Reusable Connector Cable)
 - 2130ac (Disposable Radiofrequency Grounding Pad)

Suggested Supplies:

- Infiltration solution (see below)
- Local (skin) anesthetic medication *per provider preference*

- Assortment of needles: 16g Nokor, 18g, 27g, 30g *per provider preference*
- Assortment of syringes: TB/1cc (leur lock), 3cc, 5cc, 10cc, 20cc *per provider preference*
- Infiltration Cannula
- Sterile gauze - 4x4's
- Alcohol prep pads and prep solutions such as Hibiclens or Betadine
- Towel clamps
- Sterile drapes and/or towels
- Skin marker
- Sterile gloves

Recommended for Infiltration

- Normal Saline 100 mL-1000 L bags
- 1% Lidocaine Plain 500mg/50cc bottles
- Epinephrine (1:1000) 1 mg/mL
- Sodium Bicarbonate 8.4% 1 mEq/mL

For nerve lesion procedures; Infusion Anesthetic Recommended: 4% Plain Lidocaine